


product information



product name: BAMBOO® MIS single use instrument set
in combination with: evospine MIS pedicle screw system
area: Thoraco-lumbar spine Th9 to S1
surgical approach: MIS percutaneous | Wiltse mini open

benefits

- ready to use, no waiting time
- always complete and new
- optimized daily workflow
- reduced risk for patients
- compact design/ easy handling
- works in combination with proven evospine screw



BAMBOO® set

SR-TB34-0100	instrument set Bamboo, Disposable Spinal Fixation sterile	1
	consists of:	quantity
SR-TB34-0100-01	k-wire	4
SR-TB34-0100-02	dilator	1
SR-TB34-0100-03	outer sleeve	4
SR-TB34-0100-04	screw driver	1
SR-TB34-0100-05	ratchet handle	1
SR-TB34-0100-06	trial rod 6.0x60mm	1
SR-TB34-0100-07	rod holder	1
SR-TB34-0100-08	multi-functional handle	1
SR-TB34-0100-09	break-off-nut holder	1
SR-TB34-0100-10	inner sleeve	2
SR-TB34-0100-11	torx screw driver	1
SR-TB34-0100-12	tab breaker/ counter torque	1

BAMBOO® extra parts- single packaged

SR-TB34-0100-02-01-S	Bamboo screw tap, 5,5mm, sterile	1
SR-TB34-0100-04-S	Bamboo screw driver - sterile	1
SR-TB34-0100-05-S	Bamboo ratchet handle - sterile	1
SR-TB34-0100-09-S	Bamboo break-of-nut-holder - sterile	1
SR-TB34-0100-10-01-S	Bamboo distractor - sterile	1
SR-TB34-0100-11-S	Bamboo torx screw driver - sterile	1

break-off-nut

SR-B216-01095	Bamboo break-off-nuts for fixation of evospine titanium rods
---------------	--

The dark grey material is made of high performance plastic (polyarylamid) Ixef GS-1022.
It is radiolucent and offers an optimized view during the surgery.

