

# Product information

## Bone cement and more



### Overview of bone cements:

Name	BonOs	Mendec Spine	Mendec Spine MV	Mendec Spine HV
viscosity	medium	low	medium	high
time*	5,5 / 7,5 / 9	5 / 3 / 9	3 / 8 / 10	1 / 14 / 4
amount	20,5 ml	24,5 ml	21,9 ml	21,9 ml
use cases	<ul style="list-style-type: none"> <li>• screw augmentation</li> <li>• kyphoplasty</li> <li>• vertebroplasty</li> <li>• vesselplasty</li> </ul>			
mixing system	bowl or mixer (not included)	shaker (included)	closed (included)	closed (included)

\*(preparation time / application time / polymerization time in minutes)

#### Mendec Spine



#### BonOs Inject



#### Mendec MV



#### Mendec HV



# Product information

## Bone cement and more



### Order information

Name	Ref.-No.	Description
<b>Bone cements</b>		
BonOs Inject	BS-01-0310-E	medium viscosity bone cement BonOs Inject 1 x 24 g, single
Mendec Spine	TC-1230	predosed, low viscosity, radiopaque acrylic resin with mixing system (shaking flask)
Mendec Spine Kit	TC-13C2000	radiopaque acrylic resin with applicator and tube
Mendec Spine MV	TC-13C2050	medium viscosity, acrylic resin, enclosed mixing system with application tube
Mendec Spine HV	TC-13C2040	high viscosity, acrylic resin, enclosed mixing system with application adapter
<b>Application adapters</b>		
B.Braun Combifix Adapter	WE-BR5206642	B.Braun Combifix male/male adapter, LuerLock
Tube	TC-ASB0250	male/male-tube-adapter 300mm length, Ø(outer) 2,8mm
V-Fill	TC-ASC0000	Mendec V Fill divider incl. 5 LuerLock-sealing caps, 1 LuerLock-connector and 6 caps
Mendec Aqua	TC-ASC0020	Mendec Aqua cement application system, 1 cartridge contained
Cement cartridge for Mendec Aqua	TC-ASC0030	Cement cartridge for Mendec Aqua system
<b>Mixing systems</b>		
Mini Mix	SMLV	Cement mixing system
Mixing bowl	DH-0708211000	Cement mixing bowl (400 ml) with spatula
<b>Syringes</b>		
Inject single use syringe	WE-BR04606701V	Inject single use syringe 2 ml, LuerLock
Inject single use syringe	WE-BR04606027V	Inject single use syringe 2 ml
LuerLock syringe	WE-BD 301229	LuerLock syringe, 30 ml

\* Additional items for kyphoplasty or vertebroplasty can be found in separate product information sheets „KMC-MAXX“ und „Needles“