

# Product information




Product name: cervical plate  
Area impacted: cervical spine  
Surgical approach: cervical fusion

## Features/Design

1. Polyaxial screwing
2. Backout-lock via rotation of locking screw (middle)
3. 25 to 65 mm | • additional sizes on request


## Dimensions



1-level



2-level



polyaxial | Ø 4,0 mm | screw head sinks into plate

Length is measured from the tip to the beginning of the head.  
The screw head has a length of additional 2 mm per screw.

## Implant information

### Screw

| Ref.-No.        | Description  |
|-----------------|--|
| SR-B11J-24012-S | cervical screw, variable with centred tip Ø 4 x 12 mm, light blue, sterile |
| SR-B11J-24014-S | cervical screw, variable with centred tip Ø 4 x 14 mm, pink, sterile       |
| SR-B11J-24016-S | cervical screw, variable with centred tip Ø 4 x 16 mm, silver, sterile     |
| SR-B11J-24018-S | cervical screw, variable with centred tip Ø 4 x 18 mm, dark blue, sterile  |
| SR-B11J-24515-S | cervical screw, variable with centred tip Ø 4,5 x 15 mm, gold, sterile     |

● only on request

### Plate

| Ref.-No.        | Description                    | Area    |
|-----------------|--------------------------------|---------|
| SR-B112-20025-S | cervical plate 25 mm, sterile  | 1 Level |
| SR-B112-20027-S | cervical plate 27 mm, sterile  | 1 Level |
| SR-B112-20030-S | cervical plate 30 mm, sterile  | 1 Level |
| SR-B112-20037-S | cervical plate 37 mm, sterile  | 2 Level |
| SR-B112-20040-S | cervical plate 40 mm, sterile  | 2 Level |
| SR-B112-20042-S | cervical plate 42 mm, sterile  | 2 Level |
| SR-B112-20045-S | cervical plate 45 mm, sterile  | 2 Level |
| SR-B112-20047-S | cervical plate 47 mm, sterile  | 2 Level |
| SR-B112-20055-S | cervical plate 55 mm, sterile  | 3 Level |
| SR-B112-20057-S | cervical plate 57 mm, sterile  | 3 Level |
| SR-B112-20060-S | cervical plate 60 mm, sterile  | 3 Level |
| SR-B112-20065-S | cervical plate 65 mm, sterile  | 3 Level |
| SR-B112-20067-S | cervical plate 67 mm, sterile  | 3 Level |
| SR-B112-20070-S | cervical plate 70 mm, sterile  | 3 Level |
| SR-B112-20072-S | cervical plate 72 mm, sterile  | 4 Level |
| SR-B112-20075-S | cervical plate 75 mm, sterile  | 4 Level |
| SR-B112-20077-S | cervical plate 77 mm, sterile  | 4 Level |
| SR-B112-20080-S | cervical plate 80 mm, sterile  | 4 Level |
| SR-B112-20082-S | cervical plate 82 mm, sterile  | 4 Level |
| SR-B112-20085-S | cervical plate 85 mm, sterile  | 4 Level |
| SR-B112-20090-S | cervical plate 90 mm, sterile  | 5 Level |
| SR-B112-20095-S | cervical plate 95 mm, sterile  | 5 Level |
| SR-B112-20100-S | cervical plate 100 mm, sterile | 5 Level |
| SR-B112-20105-S | cervical plate 105 mm, sterile | 5 Level |
| SR-B112-20110-S | cervical plate 110 mm, sterile | 5 Level |

● only on request